



3Y SERIES MICROBALANCES - THE NOVELTIES

- ✓ Higher resolution – up to 600 million intervals
- ✓ More precise temperature measurement
- ✓ Brand new signal filtering algorithm, enabling selective tuning to actual interfering frequency.
- ✓ Modernized mechanics design - Susceptibility to air drafts reduced six times
- ✓ Cooperation with THB module
- ✓ Brand new, faster terminal comprising: audio module (audio readout of the weighing result), video module (tutorial videos playback), WiFi interface and possibility of cooperation with applications based on ANDROID system.


New generation of microbalances MYA 3Y are intended to measure mass with the highest accuracy. They feature 5,7" LCD colour touch panel which provides new possibilities of balance operation and presenting measurement results. Measurement reliability and accuracy are maintained by system of automatic internal adjustment/calibration.

Level control in the MYA 3Y series is based on LevelSENSING system, RADWAG patented solution, which uses a system of an electronic level. New function is online monitoring of ambient conditions through built-in ambient conditions module with visualization on balance's display. Personalization of balance settings is carried out in extended user profiles and multilevel system of access levels to balance's menu. Control over opening and closing of the weighing chamber's door is carried out through programmable IR proximity sensors, installed on the terminal's overlay.

Differential weighing mode aids mass control of the same sample subjected to differed processes over time. It is particularly useful in pharmacy, environmental protection, petroleum chemistry, etc. Pipettes calibration mode is carried out using an adapter, which is an ergonomic tool aiding calibration and checking of piston pipettes using gravimetric measuring method. Extended databases enable storing all carried out measurements, with option of printing and exporting them. Standard and user defined printouts allow for maintaining documentation complying with GLP/GMP requirements practically in any application.

- Filling
- Checkweighing
- Percent setup
- Statistics
- Infrared sensors
- Air buoyancy compensation
- Autotest
- Labelling
- GLP procedures
- Formulas

Dimensions:



Technical data:

	MYA 2.3Y	MYA 0,8/3.3Y	MYA 5.3Y	MYA 11.3Y	MYA 21.3Y
Max load	2 g	0,8 g / 3 g	5 g	11 g	21 g
Readability	1 µg	1 µg / 10 µg	1 µg	1 µg	1 µg
Repeatability *	1 µg (2g)	1 µg (0,8g) 5 µg (0,8g÷3g)	1 µg (2g) 1,6 µg (2g÷5g)	1,5 µg (to 0,2g) 2,0 µg (0,2g÷5g) 2,5 µg (5g÷11g)	1,5 µg (to 0,2g) 2,0 µg (0,2g÷5g) 2,5 µg (5g÷11g) 3,0 µg (11g÷21g)
Linearity	±3 µg	±3 µg / ±4 µg	±5 µg	±6 µg	±7 µg
Eccentric load deviation	3 µg	3 µg / 4 µg	5 µg	6 µg	7 µg
Sensitivity offset	$1,5 \times 10^{-6} \times Rt$	$1,5 \times 10^{-6} \times Rt$	$1,5 \times 10^{-6} \times Rt$	$3 \times 10^{-6} \times Rt$	$4 \times 10^{-6} \times Rt$
Sensitivity temperature drift	$1 \times 10^{-6} / ^\circ C \times Rt$	$1 \times 10^{-6} / ^\circ C \times Rt$	$1 \times 10^{-6} / ^\circ C \times Rt$	$1 \times 10^{-6} / ^\circ C \times Rt$	$1 \times 10^{-6} / ^\circ C \times Rt$
Sensitivity stability	$1 \times 10^{-6} / Rok \times Rt$	$1 \times 10^{-6} / Rok \times Rt$	$1 \times 10^{-6} / Rok \times Rt$	$1 \times 10^{-6} / Rok \times Rt$	$1 \times 10^{-6} / Rok \times Rt$
Minimum weight (USP)	2 mg	2 mg	2 mg	3,0 mg	3,0 mg
Minimum weight (U = 1%, k = 2)	0,2 mg	0,2 mg	0,2 mg	0,3 mg	0,3 mg
Pan size	∅ 16 mm	∅ 16 + 60 mm (do filtrów)	∅ 26 mm	∅ 26 mm	∅ 26 mm
Weighing chamber dimensions	∅ 90 × 90 mm				
Stabilization time	5 s				
Adjustment/Calibration	automatic (internal)				
Power supply	13,5 ÷ 16 V DC / 2,1 A				
Casing of the terminal	ABS plastic				
Display	colour 5,7"(640x480) with a resistive touch screen				
Processor	2 × 1 GHz				
Memory	RAM: 256 MB DDR2, flash: 8 GB microSD				
Interface	2×USB host, 2×RS 232, Ethernet 10/100 Mbit, WiFi 802.11 b,g,n - optional				
Audio module	YES (voice messages support)				
Video support	YES (videos and multimedia instructions)				
IN / OUT	4 in / 4 out (digital)				
Ambient conditions					
Working temperature	+10 ° ÷ +40 °C				
Change rate of working temperature	±0,3 °C/h (±1 °C/8h)				
Atmospheric humidity	40% ÷ 80%				
Change rate of atmospheric humidity	±1%/h (±4%/8h)				

Rt - net weight

* Repeatability is expressed as a standard deviation from 10 weighing cycles.

Data given in tables are values determined in typical laboratory conditions. In the actual operation conditions the values of these parameters may differ from those listed above because of the impact of ambient conditions and/or balance settings.

Additional equipment:

Antivibration table for microbalances	Antistatic ionizer DJ-03
Professional weighing table	THB 2 ambient conditions module
Impact Epson printer	Additional LCD display "WD-5"
Label printer Citizen	PC USB keyboard
Anti draft shield for microbalances	Power adapter with battery and charger ZR-02
Tare and Print foot button	Mass standard
PW-WIN computer software	Antistatic cable PA 1
RAD-KEY computer software	Bar code scanner
REC-FS computer software	Cable RS 232 (balance - Epson, Citizen printer) "P0151"



New generation of microbalances MYA 3Y is designed to meet the highest requirements for determination of mass. Measurement reliability and accuracy are maintained by system of automatic internal adjustment / calibration.

Microbalances comprise two major components (an indicator and a precise mechanical measuring system are enclosed separately). Such design eliminates the influence of heat sourcing from instrument's electronics on its mechanical components and additionally protects it from shocks and vibrations caused by users operating the instrument.

All the elements of a microbalance are made of glass and steel which eliminates the influence of electrostatics on weighing process.


3Y SERIES MICROBALANCES - THE NOVELTIES

- ✓ Higher resolution – up to 600 million intervals
- ✓ More precise temperature measurement
- ✓ Brand new signal filtering algorithm, enabling selective tuning to actual interfering frequency.
- ✓ Modernized mechanics design - Susceptibility to air drafts reduced six times
- ✓ Cooperation with THB module
- ✓ Brand new, faster terminal comprising: audio module (audio readout of the weighing result), video module (tutorial videos playback), WiFi interface and possibility of cooperation with applications based on ANDROID system.


- Filling
- Checkweighing
- Percent setup
- Statistics
- Infrared sensors
- Air buoyancy compensation
- Autotest
- Labelling
- GLP procedures
- Formulas

Dimensions:

F version



F1 version



Technical data:

	MYA 0,8/3.3Y	MYA 5.3Y.F	MYA 5.3Y.F1
	-	-	-
Max load	0,8/3 g	5 g	5 g
Readability	1/10 µg	1 µg	1 µg
Repeatability *	1 µg / 5 µg	1,6 µg (< 2 g); 2,5 µg (2 g ÷ 5 g)	1,6 µg (< 2 g); 2,5 µg (2 g ÷ 5 g)
Linearity	±3 µg / ±4 µg	±5 µg	±5 µg
Eccentric load deviation	3 µg / 4 µg	5 µg	5 µg
Sensitivity offset	$1,5 \times 10^{-6} \times Rt$	$1,5 \times 10^{-6} \times Rt$	$1,5 \times 10^{-6} \times Rt$
Sensitivity temperature drift	$1 \times 10^{-6} / ^\circ C \times Rt$	$1 \times 10^{-6} / ^\circ C \times Rt$	$1 \times 10^{-6} / ^\circ C \times Rt$
Sensitivity stability	$1 \times 10^{-6} / Rok \times Rt$	$1 \times 10^{-6} / Rok \times Rt$	$1 \times 10^{-6} / Rok \times Rt$
Minimum weight (USP)	2 mg	3,2 mg	3,2 mg
Minimum weight (U = 1%, k = 2)	0,2 mg	0,32 mg	0,32 mg
Pan size	∅16 + ∅60 mm (weighing pan for filters)	∅ 100 mm + ∅ 26 mm	∅ 160 mm + ∅ 26 mm
Weighing chamber dimensions	∅ 90 × 90 mm	∅ 118 × 35 mm	∅ 168 × 35 mm
Stabilization time		5 s	
Adjustment/Calibration		automatic (internal)	
Power supply		13,5 ÷ 16 V DC / 2,1 A	
Casing of the terminal		ABS plastic	
Display		colour 5,7"(640x480) with a resistive touch screen	
Processor		2 × 1 GHz	
Memory		RAM: 256 MB DDR2, flash: 8 GB microSD	
Interface		2×USB host, 2×RS 232, Ethernet 10/100 Mbit, WiFi 802.11 b,g,n - optional	
Audio module		YES (voice messages support)	
Video support		YES (videos and multimedia instructions)	
IN / OUT		4 in / 4 out (digital)	
Ambient conditions			
Working temperature		+10 ° ÷ +40 °C	
Change rate of working temperature		±0,3 °C/h (±1 °C/8h)	
Atmospheric humidity **		40% ÷ 80%	
Change rate of atmospheric humidity		±1%/h (±4%/8h)	

Rt - net weight

* - repeatability expressed as standard deviation from 10 weighing cycles

** - Non-condensing conditions

Additional equipment:

Antivibration table for microbalances	Antistatic ionizer DJ-03
Professional weighing table	THB 2 ambient conditions module
Impact Epson printer	Additional LCD display "WD-5"
Label printer Citizen	PC USB keyboard
Anti draft shield for microbalances	Power adapter with battery and charger ZR-02
Tare and Print foot button	Mass standard
PW-WIN computer software	Antistatic cable PA 1
RAD-KEY computer software	Bar code scanner
REC-FS computer software	Cable RS 232 (balance - Epson, Citizen printer) "P0151"












MYA 3Y.P Microbalances for calibration of pipettes



release date 12-08-2014



A - evaporation ring
B - calibration vessel

-  Filling
-  Checkweighing
-  Percent setup
-  Statistics
-  Infrared sensors
-  Air buoyancy compensation
-  Autotest
-  Labelling
-  GLP procedures
-  Formulas
-  Pipette calibration


3Y SERIES MICROBALANCES - THE NOVELTIES

- ✓ Higher resolution – up to 600 million intervals
- ✓ More precise temperature measurement
- ✓ Brand new signal filtering algorithm, enabling selective tuning to actual interfering frequency.
- ✓ Modernized mechanics design - Susceptibility to air drafts reduced six times
- ✓ Cooperation with THB module
- ✓ Brand new, faster terminal comprising: audio module (audio readout of the weighing result), video module (tutorial videos playback), WiFi interface and possibility of cooperation with applications based on ANDROID system.

New generation of microbalances MYA 3Y are designed to meet the highest requirements of mass measurements. Balance's weighing chamber is adjusted to calibration of piston pipettes. The non-central location of the opening in chamber's top cover facilitates pipette insertion. Measurement reliability and accuracy are maintained by system of automatic internal adjustment / calibration. Microbalances consist of two major parts (an indicator and a precise mechanical measurement system in a separate enclosure). This solution eliminates the temperature influence and separates from shocks and vibrations caused by users operating software. All the elements of the balance are made of glass and steel which eliminates the influence of electrostatics on the weighing process.

Additional adapter for pipettes calibration is a standard equipment of the balance.

Dimensions:



Technical data:

MYA 21.3Y.P

M

Max load	21 g
Readability	1 µg
Repeatability *	1,5 µg (to 0,2g) 2,0 µg (0,2g÷5g) 2,5 µg (5g÷11g) 3,0 µg (11g÷21g)
Linearity	±7 µg
Eccentric load deviation	7 µg
Sensitivity offset	$4 \times 10^{-6} \times R_t$
Sensitivity temperature drift	$1 \times 10^{-6} / ^\circ\text{C} \times R_t$
Sensitivity time drift	$1 \times 10^{-6} / \text{Rok} \times R_t$
Minimum weight (USP)	3,0 mg
Minimum weight (U = 1%, k = 2)	0,3 mg
Pan size	∅ 26 mm
Weighing chamber dimensions	∅ 90 × 90 mm
Stabilization time	5 s
Adjustment / Calibration	automatic (internal)
Working temperature	+10 ° ÷ +40 °C
Relative air humidity **	40% ÷ 80%
Power supply	13,5 ÷ 16 V DC / 2,1 A
Casing of the terminal	ABS plastic
Display	colour 5,7"(640x480) with a resistive touch screen
Processor	2 × 1 GHz
Memory	RAM: 256 MB DDR2, flash: 8 GB microSD
Interface	2×USB host, 2×RS 232, Ethernet 10/100 Mbit, WiFi 802.11 b,g,n - optional
Audio module	YES (voice messages support)
Video support	YES (videos and multimedia instructions)
IN / OUT	4 in / 4 out (digital)

Rt - net weight

* Repeatability is expressed as a standard deviation from 10 weighing cycles

** Non-condensing conditions

Additional equipment:

Antivibration table for microbalance	Antistatic ionizer DJ-03
Profesional weighing table	Ambient conditions module
Impact Epson printer	Additional LCD display "WD-5"
Label printer Citizen	PC keyboard
Anti draft shield for microbalances	Power adapter with battery and charger ZR-02
Air density determination kit	Mass standard
Tare and "Print" foot button	Antistatic cable
PW-WIN computer software	Bar code scanner
RAD-KEY computer software	Cable RS 232 (scale - computer) "P0108"
REC-FS computer software	Cable RS 232 (scale, Epson, Citizen printer) "P0151"
Pipettes computer software	













New generation of ultra-microbalances UYA 3Y is designed to meet the high requirements for mass determination with highest accuracy of 0,1 µg and 2 g capacity.

Measurement reliability and accuracy is ensured by internal adjustment/ calibration system.

Microbalances comprise two major components (an indicator and a precise mechanical measuring system are enclosed separately). Such design eliminates the influence of heat sourcing from instrument's electronics on its mechanical components and eliminates vibration transfer.


All the elements of a microbalance are made of glass and steel which eliminates the influence of electrostatics on weighing process.

-  Filling
-  Checkweighing
-  Percent setup
-  Statistics
-  Infrared sensors
-  Air buoyancy compensation
-  Autotest
-  Labelling
-  GLP procedures
-  Formulas

3Y SERIES MICROBALANCES - THE NOVELTIES

- ✓ Higher resolution – up to 600 million intervals
- ✓ More precise temperature measurement
- ✓ Brand new signal filtering algorithm, enabling selective tuning to actual interfering frequency.
- ✓ Modernized mechanics design - Susceptibility to air drafts reduced six times
- ✓ Cooperation with THB module
- ✓ Brand new, faster terminal comprising: audio module (audio readout of the weighing result), video module (tutorial videos playback), WiFi interface and possibility of cooperation with applications based on ANDROID system.

Dimensions:



Technical data:

	UYA 2.3Y
	-
Max load	2 g
Min load	0,01 mg
Readability	0,1 µg
Repeatability	0,4 µg (0,2g÷1g) 0,6 µg (1g÷2g)
Linearity	±1,5 µg
Eccentric load deviation	1,5 µg
Sensitivity offset	$1,5 \times 10^{-6} \times R_t$
Sensitivity temperature drift	$1 \times 10^{-6} / ^\circ\text{C} \times R_t$
Sensitivitytime drift	$1 \times 10^{-6} / R_{ok} \times R_t$
Minimum weight (USP)	0,9 mg
Minimum weight (U = 1%, k = 2)	0,08 mg
Pan size	ø 16 mm
Weighing chamber dimensions	ø 90 × 90 mm
Stabilization time	10-20 s
Adjustment / Calibration	automatic (internal)
Power supply	13,5 ÷ 16 V DC / 2,1 A
Casing of the terminal	ABS plastic
Display	colour 5,7"(640x480) with a resistive touch screen
Processor	2 × 1 GHz
Memory	RAM: 256 MB DDR2, flash: 8 GB microSD
Interface	2×USB host, 2×RS 232, Ethernet 10/100 Mbit, WiFi 802.11 b,g,n - optional
Audio module	YES (voice messages support)
Video support	YES (videos and multimedia instructions)
IN / OUT	4 in / 4 out (digital)
Ambient conditions:	
Working temperature	+18 ° ÷ +30 °C
Change rate of working temperatures	±0,3°C/1h (±0,5°C / 12h)
Relative air humidity **	40% ÷ 80%
Change rate of atmospheric humidity	40% ÷ 60% (±5% / 4h)

Rt - net weight

** - Non-condensing conditions

Additional equipment:

Anti vibration table for microbalance	Ambient conditions module
Profesional weighing table	Additional LCD display "WD-5"
Impact Epson printer	PC keyboard
Label printer Citizen	Power adapter with battery and charger ZR-02
Anti draft shield for microbalances	Mass standard
"Tare" and "Print" foot button	Antistatic cable
"PW-WIN" computer software	Bar code scanner
"RAD-KEY" computer software	Cable RS 232 (balance - computer) "P0108"
"REC-FS" computer software	Cable RS 232 (balance - Epson, Citizen printer) "P0151"
Antistatic ionizer DJ-03	