

ESCO

WORLD CLASS. WORLDWIDE.



CelCulture® CO₂ Incubators

CelCulture®

CO₂ Incubators
Cradle for Beautiful Cells



Designed in the USA




WELCOME TO ESCO

Esco's Vision is to provide enabling technologies for scientific discoveries to make human lives healthier and safer.

- A leader in the development of controlled environment, laboratory and pharmaceutical equipment solutions.
- A world leader in biological safety cabinets.
- Esco has established offices in 13 countries such as Bahrain, China, India, Japan, Korea, Malaysia, Philippines, Singapore, UK, US, Vietnam, South Africa and Indonesia and is continually expanding.
- North American facilities in Pennsylvania; sales, service, logistics for US & Canada.
- Group total of more than 600 employees.
- Distributors in more than 100 countries.
- Products independently tested to international standards.
- Large R&D investments, world leading technologies.
- State-of-the-art production; vertically integrated manufacturing floor space.
- Worldwide service played out over a geographic expanse so broad that the sun never sets on what we do.



GLOBAL NETWORK



PRODUCTS AND APPLICATION

Esco Life Science Tools

Laboratory Equipment

Biosafety and Laminar Flow

Class II Type A2 Biological Safety Cabinets
 Class II Type B2 Biological Safety Cabinets
 Class III Biological Safety Cabinets
 Horizontal Laminar Flow Clean Benches
 Vertical Laminar Flow Clean Benches
 Laboratory Animal Research Workstations
 PCR Cabinets

Fume Hoods

Laboratory Fume Hoods
 Ductless Fume Hoods
 Fume Scrubbers
 Exhaust Blowers
 Fume Hood Airflow Monitors

PCR

PCR Thermal Cyclers
Not Available in North America

Incubators and Ovens

Forced Convection Laboratory Ovens
 Forced Convection Laboratory Incubators
 Refrigerated Incubators
 CO₂ Incubators
 Remote Monitoring, Data Logging, and Programming Software

Cold Storage

Ultra-low Temperature Freezers

Medical Equipment

Assisted Reproductive Technology

ART Workstations
 CO₂ Incubators with Suppressed O₂
 Multi-room Incubators

Pharmaceutical Equipment

Containment / Compounding Pharmacy

Downflow Booths
 Powder Weighing Balance Enclosures
 Pharmacy Isolators
 Cytotoxic Safety Cabinets
 Soft Capsule
 Air Showers
 Straddle Units
 Garment Storage Cabinets
 Pass Boxes
 Transfer Hatches



CelCulture® CO₂ Incubators

INTRODUCTION

CO₂ incubators are widely used in scientific research to grow and maintain cell cultures. Typical fields of application include tissue engineering, *in vitro* fertilization, neuroscience, cancer research and other mammalian cell research.

Sleek, reliable and intuitive, Esco CelCulture CO₂ incubators provide all-rounded sample protection that brings your scientific dreams one step closer to reality.

KEY FEATURES

CelCulture® CO₂ INCUBATORS *Cradle for Beautiful Cells*



CelCulture® CO₂ Incubators
are available in 3 sizes, 50L, 170L, 240L.

ULPA FILTER

- 99.999% efficient, superior to conventional HEPA filters
- Filters air continuously
- Chamber returns to ISO Class 5 cleanliness in 13 minutes upon door closing to prevent contamination



SHELVING

- Perforated shelving to improve uniformity
- Anti-tip
- Stainless steel
- Built-in grip
- Dismantles without tools for easy cleaning

DIRECT HEAT & AIR JACKET

- Fast and uniform heating
- Rapid temperature recovery without overshoot
- Air jacket improves chamber stability



DUCT WORK

- Directs air flow for rapid recovery and excellent uniformity
- Easily removed for cleaning



WATER PAN

- Precisely heated by base heater to provide high humidity
- Gentle airflow over water surface accelerates humidity recovery



ROUNDED CORNERS

- Seamless design
- Facilitates cleaning

O₂ SENSOR

for suppressed O₂ model

- Long life
- Stable output signal
- No influence from CO₂

DOOR SWITCH

Automatically turns off the blower and gas supply when the door is opened

TOP COVER

Provides quick access to electrical panel components

CO₂ SENSOR

- Choice of TC or IR
- Single-beam, dual-wavelength IR sensor is drift-free
- Auto-zeroing

SMARTSENSE™ MICROPROCESSOR INTERFACE

Intuitive, fully equipped control and monitoring system

BLOWER

Gentle airflow in chamber improves recovery and uniformity

OUTER DOOR

- Reversible
- Heated to prevent condensation

SAMPLE PORT

Allows direct measurement of chamber atmosphere such as temperature and CO₂ concentration

GLASS DOOR

For observing sample cells inside the chamber during operation

DOOR LATCH

To lock / unlock the glass door

LEVELING FEET

Easily adjustable


QUALITY ESCO CONSTRUCTION

- Electrogalvanized steel with white oven-baked epoxy-polyester antimicrobial powder-coated finish.
- External surfaces are powder coated with Esco **ISOCIDE™** to eliminate 99.9% of surface bacteria within 24 hours of exposure.
- Ensures a healthier, safer and cleaner lab environment.


Note: For 50 L Model, no top plenum, bottom plenum, blower & ULPA Filter.

VIVOCELL™ PRECISE PARAMETER CONTROL

BEST UNIFORMITY AND CONTROL AMONG COMPETITION




Different lines represent different sensor positions inside the chamber. Esco CelCulture has uniformity variance of less than ± 0.2 °C which means all the samples are evenly heated.*



Minimal fluctuation (± 0.1 °C) ensures temperature stability.*



FAST CO₂, TEMPERATURE AND HUMIDITY RECOVERY WITHOUT OVERTHOOT



Precisely tuned sensor and software result in fast recovery of CO₂ without overshoot. This ensures uniform CO₂ levels even with frequent incubator door openings.* Similarly, temperature and humidity recoveries are twice as fast as conventional incubators.

- Company A's model: overshoot.
- Company B's model: slow recovery.
- Esco CelCulture: fast recovery, no overshoot.



DIRECT HEAT AND AIR JACKET



- Direct heating enables rapid temperature recovery while air jacket provides isolation against ambient temperature fluctuations.
- Precise heating in the chamber is achieved by using 8 heaters (3 zones). The 3 zones are intelligently controlled by the microprocessor for best temperature uniformity and minimal fluctuation.

- The main heater provides precise temperature control.
- The bottom heater warms the water pan and controls humidity.
- The outer door heater prevents condensation on glass door and facilitates temperature recovery.

VENTIFLOW™ FORCED CONVECTION



- No disturbance to cell culture.
- Blower automatically stops when door is opened to minimize mixing of chamber and room air.
- Accelerates recovery of chamber air to ISO Class 5 Cleanliness after door closing to prevent contamination.
- Improves CO₂, humidity and temperature uniformity.
- Filtered air circulates across water pan to accelerate humidifying process.

ORDERING INFORMATION

TC SENSOR MODEL WITH STAINLESS STEEL CHAMBER

MODELS	DESCRIPTION
CCL-050A-8	CelCulture® Incubator, 50L, TC Sensor, CO ₂ Control, Moist Heat Decon, 230VAC, 50/60HZ
CCL-050A-9	CelCulture® Incubator, 50L, TC Sensor, CO ₂ Control, Moist Heat Decon, 115VAC, 50/60HZ
CCL-170A-8	CelCulture® Incubator, 170L, TC Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 230VAC, 50/60HZ
CCL-170A-8-NF	CelCulture® Incubator, 170L, TC Sensor, CO ₂ Control, Moist Heat Decon, 230VAC, 50/60HZ, (No ULPA Filter)
CCL-170A-9	CelCulture® Incubator, 170L, TC Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 115VAC, 50/60HZ
CCL-170A-9-NF	CelCulture® Incubator, 170L, TC Sensor, CO ₂ Control, Moist Heat Decon, 115VAC, 50/60HZ, (No ULPA Filter)
CCL-240A-8	CelCulture® Incubator, 240L, TC Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 230VAC, 50/60HZ
CCL-240A-8-NF	CelCulture® Incubator, 240L, TC Sensor, CO ₂ Control, Moist Heat Decon, 230VAC, 50/60HZ, (No ULPA Filter)
CCL-240A-9	CelCulture® Incubator, 240L, TC Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 115VAC, 50/60HZ
CCL-240A-9-NF	CelCulture® Incubator, 240L, TC Sensor, CO ₂ Control, Moist Heat Decon, 115VAC, 50/60HZ, (No ULPA Filter)

IR SENSOR MODEL WITH STAINLESS STEEL CHAMBER

MODELS	DESCRIPTION
CCL-050B-8	CelCulture® Incubator, 50L, IR Sensor, CO ₂ Control, Moist Heat Decon, 230VAC, 50/60HZ
CCL-050B-9	CelCulture® Incubator, 50L, IR Sensor, CO ₂ Control, Moist Heat Decon, 115VAC, 50/60HZ
CCL-170B-8	CelCulture® Incubator, 170L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 230VAC, 50/60HZ
CCL-170B-8-NF	CelCulture® Incubator, 170L, IR Sensor, CO ₂ Control, Moist Heat Decon, 230VAC, 50/60HZ, (No ULPA Filter)
CCL-170B-9	CelCulture® Incubator, 170L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 115VAC, 50/60HZ
CCL-170B-9-NF	CelCulture® Incubator, 170L, IR Sensor, CO ₂ Control, Moist Heat Decon, 115VAC, 50/60HZ, (No ULPA Filter)
CCL-240B-8	CelCulture® Incubator, 240L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 230VAC, 50/60HZ
CCL-240B-8-NF	CelCulture® Incubator, 240L, IR Sensor, CO ₂ Control, Moist Heat Decon, 230VAC, 50/60HZ, (No ULPA Filter)
CCL-240B-9	CelCulture® Incubator, 240L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 115VAC, 50/60HZ
CCL-240B-9-NF	CelCulture® Incubator, 240L, IR Sensor, CO ₂ Control, Moist Heat Decon, 115VAC, 50/60HZ, (No ULPA Filter)

SURPRESSED O₂ MODEL WITH STAINLESS STEEL CHAMBER

MODELS	DESCRIPTION
CCL-050T-8	CelCulture® Incubator, 50L, IR Sensor, CO ₂ Control, O ₂ Control, Moist Heat Decon, 230VAC 50/60HZ
CCL-050T-9	CelCulture® Incubator, 50L, IR Sensor, CO ₂ Control, O ₂ Control, Moist Heat Decon, 115VAC 50/60HZ
CCL-170T-8	CelCulture® Incubator, 170L, IR Sensor, CO ₂ Control, O ₂ Control, ULPA, Moist Heat Decon, 230VAC 50/60HZ
CCL-170T-8-NF	CelCulture® Incubator, 170L, IR Sensor, CO ₂ Control, O ₂ Control, Moist Heat Decon, 230VAC 50/60HZ, (No ULPA Filter)
CCL-170T-9	CelCulture® Incubator, 170L, IR Sensor, CO ₂ Control, O ₂ Control, ULPA, Moist Heat Decon, 115VAC 50/60HZ
CCL-170T-9-NF	CelCulture® Incubator, 170L, IR Sensor, CO ₂ Control, O ₂ Control, Moist Heat Decon, 115VAC 50/60HZ, (No ULPA Filter)
CCL-240T-8	CelCulture® Incubator, 240L, IR Sensor, CO ₂ Control, O ₂ Control, ULPA, Moist Heat Decon, 230VAC 50/60HZ
CCL-240T-8-NF	CelCulture® Incubator, 240L, IR Sensor, CO ₂ Control, O ₂ Control, Moist Heat Decon, 230VAC 50/60HZ, (No ULPA Filter)
CCL-240T-9	CelCulture® Incubator, 240L, IR Sensor, CO ₂ Control, O ₂ Control, ULPA, Moist Heat Decon, 115VAC 50/60HZ
CCL-240T-9-NF	CelCulture® Incubator, 240L, IR Sensor, CO ₂ Control, O ₂ Control, Moist Heat Decon, 115VAC 50/60HZ, (No ULPA Filter)



CelMate®

CO₂ Incubators

INTRODUCTION

Esco now offers the new CelMate, 170-liter and 240-liter, entry-level cell culture CO₂ incubator with superb contamination control.

This is for customers looking for a CO₂ incubator that can provide the best protection for their cell cultures but with a limited budget.



CONTAMINATION CONTROL SYSTEMS:

ULPA Filtration System , Validated Swiftcon Overnight Decontamination Cycle (20 hours), Filtered Gas Injection Lines , and Exterior is coated with Isocide.

OPTIONS AND ACCESSORIES:

All options and accessories for CelCulture CO₂ incubators are also appropriate to use on CelMate CO₂ incubator.

CelMate® CO₂ Incubators are available in 170L and 240L.

CelMate is using TC CO₂ sensor and no decontamination pump.

ORDERING INFORMATION

MODELS	DESCRIPTION
CLM-170-A-8	CelMate® Incubator, 170L, TC Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 230VAC, 50/60HZ (Without Decon Pump)
CLM-170-A-9	CelMate® Incubator, 170L, TC Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 115VAC, 50/60HZ (Without Decon Pump)
CLM-240-A-8	CelMate® Incubator, 240L, TC Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 230VAC, 50/60HZ (Without Decon Pump)
CLM-240-A-9	CelMate® Incubator, 240L, TC Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 115VAC, 50/60HZ (Without Decon Pump)



CelCulture®

CO₂ Incubators with Copper Interior Chamber

Copper has been known for millennia to have anti-microbial properties. 100% pure solid copper interiors offer additional protection for your precious samples.

ACCESSORIES

COA-2026-F Extra Shelf (50L, Solid Copper)


Each CelCulture CO₂ Incubator comes standard with 3 shelves for 50L and it can accommodate up to a maximum of 4 shelves for 50L. Extra shelves are available.

COA-2027-F Extra Shelf (170L, Solid Copper)

Each CelCulture CO₂ Incubator comes standard with 4 shelves for 170L and it can accommodate up to a maximum of 7 shelves for 170L. Extra shelves are available.

COA-2028-F Extra Shelf (240L, Solid Copper)

Each CelCulture CO₂ Incubator comes standard with 4 shelves for 240L and it can accommodate up to a maximum of 7 shelves for 240L. Extra shelves are available.



Other options and accessories for CelCulture CO₂ incubators except for the shelves are also appropriate to use on CelCulture CO₂ incubator with Copper Interior Option.

ORDERING INFORMATION

IR SENSOR MODEL WITH 100% COPPER CHAMBER


MODELS	DESCRIPTION
CCL-050B-8-Cu	CelCulture® Incubator, 50L, IR Sensor, CO ₂ Control, Moist Heat Decon, 230VAC, 50/60HZ
CCL-050B-9-Cu	CelCulture® Incubator, 50L, IR Sensor, CO ₂ Control, Moist Heat Decon, 115VAC, 50/60HZ
CCL-170B-8-Cu	CelCulture® Incubator, 170L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 230VAC, 50/60HZ
CCL-170B-9-Cu	CelCulture® Incubator, 170L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 115VAC, 50/60HZ
CCL-240B-8-Cu	CelCulture® Incubator, 240L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 230VAC, 50/60HZ
CCL-240B-9-Cu	CelCulture® Incubator, 240L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, 115VAC, 50/60HZ

SURPRESSED O₂ MODEL WITH 100% COPPER CHAMBER

MODELS	DESCRIPTION
CCL-050T-8-Cu	CelCulture® Incubator, 50L, IR Sensor, CO ₂ Control, O ₂ Control, Moist Heat Decon, 230VAC 50/60HZ
CCL-050T-9-Cu	CelCulture® Incubator, 50L, IR Sensor, CO ₂ Control, O ₂ Control, Moist Heat Decon, 115VAC 50/60HZ
CCL-170T-8-Cu	CelCulture® Incubator, 170L, IR Sensor, CO ₂ Control, O ₂ Control ULPA, Moist Heat Decon, 230VAC 50/60HZ
CCL-170T-9-Cu	CelCulture® Incubator, 170L, IR Sensor, CO ₂ Control, O ₂ Control ULPA, Moist Heat Decon, 115VAC 50/60HZ
CCL-240T-8-Cu	CelCulture® Incubator, 240L, IR Sensor, CO ₂ Control, O ₂ Control ULPA, Moist Heat Decon, 230VAC 50/60HZ
CCL-240T-9-Cu	CelCulture® Incubator, 240L, IR Sensor, CO ₂ Control, O ₂ Control ULPA, Moist Heat Decon, 115VAC 50/60HZ

CONTROLLER TYPE

USER - FRIENDLY SOFTWARE INTERFACE



- 1. Start / stop decontamination cycle
- 2. **HEAT LED**
Lights when heat is applied to chamber
- 3. °C is lit when displaying the temperature
- 4. %RH is lit when displaying the humidity level
- 5. Enter menu / go back to previous menu
- 6. Scroll up / increase value
- 7. **ALARMS LED**
Will blink when errors and warnings occur
- 8. Mute alarms
- 9. **INJECT LED**
Lights when gas is injected
- 10. %O₂ is lit when displaying the O₂ concentration
- 11. %CO₂ is lit when displaying the CO₂ concentration
- 12. Confirm value / enter a menu
- 13. Scroll down / decrease value

• **Comprehensive, user-configurable alarms:**

- Temperature
- CO₂
- Humidity (if installed)
- O₂ (if installed)

• **CelAlert™ alarm system reminds user to replace parts.**




CO₂ tank depletion reminder in addition to CO₂ tank low alarm. Automatic calculation of how much CO₂ gas is left in the tank provides fail proof reminder that alerts user one week before the gas is depleted. This gives user some buffer time to place order for new tanks.



ULPA reminder will alert user to replace ULPA filter.

• **Intelligent data and event logger records all incubator parameters for on screen recall. A 16Mb built-in flash memory guarantees long term storage of data.**



• **Diagnostic interface and online quick help provide comprehensive solutions to frequently encountered problems.**



Voyager®

Remote Monitoring, Datalogging, Programming Software

Esco Voyager® is a PC-based software package developed for the remote monitoring, datalogging, and programming / device configuration of Esco thermostatic products


It is a centralized monitoring and control system for your laboratory which provides EXTRA PROTECTION FOR YOUR SAMPLES.

Voyager® interfaces with individual Esco equipment over RS485 using the EscoBUS communications protocol. Multiple equipment maybe interfaced to a single PC.

Compatible Equipment

- Lexicon® – Ultra-low Temperature Freezer
- CelCulture® – CO₂ Incubator (CCL)
- CelMate® - CO₂ Incubator (CLM)
- Isotherm® – Forced Convection Oven (OFA)
- Isotherm® – Forced Convection Incubator (IFA)
- Isotherm® – Refrigerated Incubator (IFC)

REAR PANEL



1 Power Supply Inlet

The power supply inlet connects the incubator unit to the power source.



2 Cooling Fan

The cooling fan prevents the electrical panel from overheating.



3 RS485 Communication Port

The RS485 provides serial communication port for PC. It can be daisy chained from product to product and connected to a PC.



4 Analog Port (Optional)

The analog port allows the incubator to output analog signals representing temperature, CO₂/O₂* concentration and relative humidity, depending on the options available in the incubator. This allows the Incubator to be connected to an in-house data acquisition or alarm system.



5 Alarm Contact

A set of relay contacts located on the rear of the unit is provided to monitor temperature, humidity or CO₂ alarms. The alarm contacts can be connected to a remote alarm system.



6 N₂ Gas Supply Inlet (for Suppressed O₂ model)

The N₂ gas supply inlet is only applicable for models with N₂* Control function. Inlet pressure requirement is 15 psi.



7 CO₂ Gas Supply Inlet

The CO₂ gas supply inlet connects the CO₂ gas supply with the Incubator unit. Inlet pressure requirement is 15 psi.



8 Gas Inline Filter

Inline filters are provided to remove any contaminants from the gas supply.



9 Access Port

Allows cables, hoses or additional sensors to be routed into the work space. A rubber stopper with controlled leak is installed as standard configuration and is part of standard accessories.

* O₂ and N₂ functions are applicable only to models with Suppressed O₂.

CelCulture CO₂ INCUBATOR SENSORS



IR SENSOR

An IR sensor is a versatile instrument for measuring CO₂ level inside the incubator. The CARBOCAP® sensor is silicon based and its operation is based on the NDIR Single-Beam Dual-Wavelength principle. IR based sensors are not affected by water vapor, dust or most chemicals. The single-beam dual-wavelength technology (one reference and one measurement) ensures a drift-free sensor that does not require calibration by the user.

Operating principle

The light source is positioned to shine at the IR detector so that the light travels a fixed distance to the detector, where the intensity of the light is measured. A Fabry-Perot Interferometer (FPI) is positioned just in front of the IR detector. The FPI is a tunable filter which allows only certain wavelengths of light to pass through to the detector.

Carbon dioxide absorbs certain wavelengths of light and not others, so the FPI is designed to pass light at a CO₂ absorption wavelength (4.26 μm) and a nearby, non-absorbing wavelength.

When the sensor is operating, the FPI is regularly tuned back and forth between the two wavelengths. At the CO₂ absorption wavelength, the intensity of detected light is reduced in proportion to the concentration of CO₂ in the optical path. The light intensity measured at the non-absorbing wavelength serves as a baseline for comparison.

Operating Conditions:

%CO₂ detection range: 0 to 20% CO₂ Concentration

%RH operating range: Not affected by Humidity

Temperature range: -20°C to +60°C



TC CO₂ SENSOR

Esco TC CO₂ sensor's operating principle relies on a resistor as a heater and two thermocouples as a sensing element for the CO₂ gas. Accurate sensing is made possible by the porous cap on the eye of the sensor probe.

One of the thermocouples functions as a reference signal, while the other functions as the sensing signal. An amplifier will feed the data variance between the two thermocouples to an electronic control system.

Operating Conditions:

%CO₂ detection range: 0 to 20% CO₂ Concentration

%RH operating range: 40% to 98% Relative Humidity

Temperature range: +25°C to +100°C



O₂ SENSOR

Figaro's O₂ sensor is a unique galvanic cell type oxygen sensor. Its most notable features are long life expectancy, excellent chemical durability, and it is not influenced by CO₂. The O₂ sensor is ideal to meet the ever-increasing demand for oxygen monitoring in various fields such as combustion gas monitoring, the biochemical field, medical applications, domestic combustion appliances, etc.

Operating Conditions:

%O₂ detection range: 1 to 20.7% O₂ Concentration

%RH operating range: 10% to 90% Relative Humidity

Temperature range: 5°C to +40°C

TESTING & CERTIFICATION



For IVF applications, EscO CelCulture CO₂ incubators are certified EMBRYO-SAFE.

Rigorously tested with the Mouse Embryo Assay (MEA), the CelCulture remarkably has 100% embryo survival. The Mouse Embryo Assay (MEA) is the de facto standard test done to demonstrate that a procedure or an article of equipment is safe to use for manipulating human embryos (e.g., *in vitro* fertilization or IVF).



The EscO CelCulture CO₂ incubators is listed by Underwriters Laboratory (UL)*, to meet the requirements of both the U.S. and Canada standards for electrical/mechanical integrity.

*applicable for 170L



HPA Validated Decontamination Cycle

The EscO CelCulture CO₂ Incubator 90°C decontamination cycle has been evaluated and shown to be an effective method for deactivation of the normally resistant fungi and bacterial spores *Aspergillus brasiliensis* and *Bacillus atrophaeus*, and the vegetative cells *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Enterobacter faecalis* and *Escherichia Coli*.

OPTIONS AND ACCESSORIES



COA-1001 / COA-1001-F Humidity Display

This option allows the Incubator to monitor the relative humidity inside the chamber. The probe for the sensor works in freezing conditions (-70°C) and also in temperatures up to +180°C. The sensor is easy to install and has excellent accuracy. The airflow in the chamber does not affect the measurement. The sensor is maintenance free. It does not need to be removed during 90°C moist heat decontamination cycle.



COA-1002 / COA-1002-F CO₂ Backup

This option allows two tanks of CO₂ to be connected to the incubator. It will automatically switch from the primary tank to the secondary tank when low gas pressure is detected on the primary tank.



COA-1005 / COA-1005-F Analog Output

A set of relay contacts are provided at the rear of the incubator that allows the incubator to output analog signals representing the temperature, CO₂ / O₂ content and relative humidity, depending on the options available in your incubator. This allows the chamber to be connected to an in-house data acquisition or alarm system. This option can also be field-installed.

The analog signal outputs can be set to operate in either voltage DC (0-5 Vdc) or current (4-20 mA) mode. The factory default setting is voltage. Switch on the board to toggle between the modes.




COA-1006 / COA-1006-F Sealed Inner Door Kit (170L)

CelCulture CO₂ Incubators can be equipped with 4 glass doors, which allows access to defined sections of the incubator without disturbing the inner atmosphere. This minimizes recovery times and contamination risks. The Sealed Inner Door is available as a factory installed option or field installed retrofit kit.



COA-1007 / COA-1007-F N₂ Back-up

This option allows two tanks of N₂ to be connected to the incubator. It will automatically switch from the primary tank to the secondary tank when low gas pressure is detected on the primary tank.



COA-2018-F (50L) / COA-2001-F (170L) / COA-2019-F (240L) Roller Base

Roller base is available with casters for mobility of your incubators and to provide protection against floor contamination.



COA-2020-F (50L) / COA-2002-F (170L) / COA-2021-F (240L) Floor Stand 200 mm (8.0") With Adjustable Feet

Floor stands are available with adjustable feet, with a nominal range of 180 mm to 250 mm (7.1" to 9.8") for comfortable access to the incubator and to avoid floor contamination.



COA-2022-F (50L) / COA-2003-F (170L) / COA-2023-F (240L) Floor Stand 700 mm (27.6") With Casters

This support stand raises the incubator to a height of 700 mm (27.6") above the floor for comfortable access. It comes with casters for mobility of your incubators.



COA-2005-F 2-Stage Gas Regulator for CO₂/N₂

CO₂ and N₂ gas input regulators reduce pressure from the tank to the incubator. It has dual pressure gauges, barbed line connection and shut-off valve. It prevents over-pressurization of the gas supply into the incubator which could cause the tubing to burst.

- CGA 320 connector (U.S. Standard)
- BP-BS341-#8-NT4 connector (British Standard)

Note: Compatible with European DIN477, French NFE29-650 and Australia AS2473

- G5/8-RH connector (China Standard)



COA-2024-F (50L) / COA-2007-F (170L) / COA-2025-F (240L) Extra Shelf (Stainless Steel) for Standard Stainless Steel Chamber

Each CelCulture CO₂ Incubator comes standard with 3 shelves for 50L / 4 shelves for 170L & 240L and it can accommodate up to a maximum of 4 shelves for 50L / 7 shelves for 170L & 240L.



COA-2008-F Stacking Kit

Stacking kit is a provision to stack one incubator on top of another incubator. Four stacking brackets are included as standard inside the Accessories Kit Box with each incubator.



COA-2010-F Electronic CO₂ Analyzer, For CO₂ / Temp Measurement

COA-2016-F Electronic CO₂ + O₂ Analyzer, For CO₂ / O₂ / Temp Measurement

COA-2017-F Electronic CO₂ + O₂ + RH Analyzer, For CO₂ / O₂ / RH / Temp Measurement

The Electronic Analyzer allows the measurement of CO₂ concentration, O₂ concentration, Relative Humidity and temperature (temperature probe already included).



COA-2012-F 6" Chart Recorder, Temp, 115/230VAC 50/60HZ

The chart recorder provides an easy-to-read graph of data vs time. It is a reliable, accurate, and stable instrument for on-the-spot written documentation of incubator chamber temperature. This model offers 6" chart of temperature data.



COA-2013-F 8" Chart Recorder, Temp/Temp, 115/230VAC 50/60HZ

The chart recorder provides an easy-to-read graph of data vs time. It is a reliable, accurate, and stable instrument for on-the-spot written documentation of incubator chamber temperature. This model offers 8" chart of temperature data and comes with 2 remote probes for dual temperature monitoring.



COA-2014-F 6" Chart Recorder, Temp/RH, 115/230VAC 50/60HZ

The chart recorder provides an easy-to-read graph of data vs time. It is a reliable, accurate, and stable instrument for on-the-spot written documentation of incubator chamber temperature. This model offers 6" chart of temperature and humidity data.



COA-2015-F Inner Door Shelving Kit (4 Sets With Total 12 Mini Shelves For One Incubator) (170L)

These mini shelves are to be used with the Sealed Inner Door Kit installed. There are 4 sets with a total of 12 mini shelves on each incubator.



5250001 Voyager Software Kit

Esco Voyager is a PC-based software package developed for the remote monitoring, datalogging and programming / device configuration of Esco controlled environment laboratory equipment. Compatible equipment includes Laboratory Ovens and Incubators, Low Temperature Incubators, CO₂ Incubators and Ultra-low Temperature Freezer.

ORDERING INFORMATION

ACCESSORIES	DESCRIPTION
COA-1001	Humidity Display, Factory Installed
COA-1001-F	Humidity Display, Field Install Kit
COA-1002	CO ₂ Backup (Tank Switcher), Factory Installed
COA-1002-F	CO ₂ Backup (Tank Switcher), Field Installed
COA-1004	Reversed Door Swing, Factory Installed
COA-1005	Analog Outputs, Factory Installed
COA-1005-F	Analog Outputs, Field Installed
COA-1006	Sealed Inner Door Kit for 170L (4 Glass Doors With Latches), Factory Installed
COA-1006-F	Sealed Inner Door Kit for 170L (4 Glass Doors With Latches), Field Installed
COA-1007	N ₂ Back-up (Tank Switcher), Factory Installed
COA-1007-F	N ₂ Back-up (Tank Switcher), Field Installed
COA-2018-F	Roller Base (50L)
COA-2001-F	Roller Base (170L)
COA-2019-F	Roller Base (240L)
COA-2020-F	Floor Stand 200 mm (8.0") With Adjustable Feet (50L)
COA-2002-F	Floor Stand 200 mm (8.0") With Adjustable Feet (170L)
COA-2021-F	Floor Stand 200 mm (8.0") With Adjustable Feet (240L)
COA-2022-F	Floor Stand 700 mm (27.6") With Casters (50L)
COA-2003-F	Floor Stand 700 mm (27.6") With Casters (170L)
COA-2023-F	Floor Stand 700 mm (27.6") With Casters (240L)
COA-2005-F	2-Stage Gas Regulator for CO ₂ /N ₂ Choose One of The Connectors Below: 1080588 - CGA 320 Connector (US Standard) 1080589 - BP-BS34-#8-NT4 Connector (British Standard) 1080590 - G5/8-RH Connector (China Standard)
COA-2024-F	Extra Shelf (50 L, Stainless Steel)
COA-2007-F	Extra Shelf (170 L, Stainless Steel)
COA-2025-F	Extra Shelf (240 L, Stainless Steel)
COA-2008-F	Stacking Kit (One Set Included With Every Unit Purchased)
COA-2010-F	Electronic CO ₂ Analyzer, For CO ₂ / Temp Measurement (With Temp. Probe)
COA-2016-F	Electronic CO ₂ + O ₂ Analyzer, For CO ₂ / O ₂ / Temp Measurement
COA-2017-F	Electronic CO ₂ + O ₂ + RH Analyzer, For CO ₂ / O ₂ / RH / Temp Measurement
COA-2011-F	IQ / OQ Documentation
COA-2012-F	6" Chart Recorder, Temp, 115/230VAC 50/60HZ
COA-2013-F	8" Chart Recorder, Temp/Temp, 115/230VAC 50/60HZ
COA-2014-F	6" Chart Recorder, Temp/RH, 115/230VAC 50/60HZ
COA-2015-F	Inner Door Shelving Kit for 170L (4 Sets With Total 12 Mini Shelves For One Incubator)
5250001	Voyager Software Kit




TECHNICAL SPECIFICATIONS

Front view




Side view

Rear view




MODEL 50L



MODEL 170L



MODEL 240L



1. Control panel*
2. On / off switch
3. Blower fan
4. ULPA filter
5. Sensors

6. Access port
7. Adjustable shelves
8. Humidity pan
9. N₂ gas supply
10. CO₂ gas supply

11. Alarm contact
12. Analog output
13. RS485
14. Cooling fan
15. Power supply inlet

*Control panel position differs for CelMate CO₂ incubator models.

GENERAL SPECIFICATIONS

CELCULTURE CO₂ INCUBATORS

	CCL-050 _ _	CCL-170 _ _ CLM-170 _ _	CCL-240 _ _ CLM-240 _ _
TEMPERATURE			
Temp. Control Method	Direct heat & air jacket using Microprocessor PID		
Temp. Range, °C	Amb. +3 to 60		
Temp. Uniformity, °C	<± 0.2*	<± 0.2*	<± 0.3*
Temp. Accuracy, °C	<± 0.1		
Recovery Time** (after 1 min. door opening, 98% from initial value)	4 mins	6 mins	6 mins
Ambient Temp. Range	18 to 34°C (64 to 93 °F)		
CO₂			
CO ₂ Control System	Microprocessor PID		
CO ₂ Range, % CO ₂	0-20		
CO ₂ Accuracy, % CO ₂	± 0.1		
CO ₂ Sensor	Infrared (IR) Sensor*** / TC Sensor		
CO ₂ Recovery Time*** (after 1 min. door opening, 98% from initial value)	Standard Unit: 8 minutes Suppressed O ₂ model: 8 minutes	Standard Unit: 4 minutes Suppressed O ₂ model: 5 minutes	Standard Unit: 5 minutes Suppressed O ₂ model: 5 minutes
O₂ SPECS (FOR SUPPRESSED O₂ MODEL)			
O ₂ Control System	Microprocessor PID		
O ₂ Range, % O ₂	1-20.7%		
O ₂ Accuracy, % O ₂	± 0.1		
O ₂ Sensor	Galvanic Cell Type		
O ₂ Recovery Time (after 1 minute door opening)	At 1.0% O ₂ by volume: 10 minutes	At 1.0% O ₂ by volume: 20 minutes	At 1.0% O ₂ by volume: 24 minutes
	At 5.0% O ₂ by volume: 6 minutes	At 5.0% O ₂ by volume: 10 minutes	At 5.0% O ₂ by volume: 12 minutes
HUMIDITY			
Humidification Method	Humidity pan		
Humidity Range, % RH	Up to 97%		
PHYSICAL CONSTRUCTION			
Interior Volume	50 L (1.8 cu.ft.)	170 L (5.7 cu.ft.)	240 L (8.5 cu.ft.)
External Dimensions (W x D x H)	500 x 500 x 655 mm (19.7" x 19.7" x 25.8")	660 x 660 x 900 mm (26.0" x 26.0" x 35.4")	750 x 665 x 900 mm (29.5" x 26.2" x 35.4")
Internal Dimensions (W x D x H)	345 x 375 x 390 mm (13.6" x 14.8" x 15.4")	505 x 530 x 635 mm (19.9" x 20.9" x 25.0")	595 x 620 x 635 mm (23.4" x 24.4" x 25.0")
Shipping Weight	70 kg (154.3 lbs)	120 kg (264.6 lbs)	155 kg (341.7 lbs)
Shipping Dimensions (W x D x H)	660 x 660 x 890 mm (26.0" x 26.0" x 35.0")	850 x 720 x 1150 mm (33.5" x 28.3" x 45.3")	860 x 830 x 1110 mm (33.9" x 32.7" x 43.7")
Number of Shelves	3	4	4
Maximum No. of Shelves	4	7	7
Shelves Area (W x D)	310 x 310 mm (12.2" x 12.2")	470 x 470 mm (18.5" x 18.5")	550 x 550 mm (21.7" x 21.7")
Max. Load per Shelf	4 kg/shelf (8.8 lbs/shelf)	11 kg/shelf (24.3 lbs/shelf)	15 kg/shelf (33.1 lbs/shelf)
Available Electrical Configuration	220 - 240 VAC, 50 / 60 Hz, 1Φ, 3.4 A		
	110 - 130 VAC, 50 / 60 Hz, 1Φ, 7.0 A		
Maximum Power Consumption	372 watts	800 watts	1100 watts
Power Consumption 37°C	37 watts	80 watts	110 watts
Interior Material	Stainless steel, type 304		
CONTAMINATION CONTROL			
Contamination Control Methods	1) Main body is electrogalvanized steel with ISOCIDE antimicrobial coating; 2) Moist 90°C OVERNIGHT decon. cycle (HPA validated); 3) 0.2 micron in-line filter for gas inputs; 4) ULPA filter****		

* Data recorded under optimum factory setting conditions

** For temperature not exceeding 37°C

*** For CO₂ not exceeding 5.2%. Recovery time with TC sensor is longer.

****Not available for 50L

ESCO

WORLD CLASS. WORLDWIDE.



CelCulture® Water-Jacketed CO₂ Incubator

CelCulture®

Water-Jacketed CO₂ Incubator
Cradle for Beautiful Cells





CelCulture®

Water-Jacketed CO₂ Incubator

INTRODUCTION

Esco CelCulture® Water-Jacketed CO₂ Incubator provides a very stable environment to grow and maintain cell cultures.

Water-Jacketed CO₂ Incubator can maintain temperature by surrounding the chamber by hot walls generated from the heated water. The heated water circulates and radiates heat around the inner chamber which maintains constant temperature.

Water-Jacketed CO₂ Incubator can hold the chamber temperature much longer when power is lost. The CO₂ incubator will also be able to recover temperature much faster after power failure as it also gets back more quickly as temperature settings change due to frequent opening of the door.

KEY FEATURES

MORE STABLE TEMPERATURE CONTROL

- Faster temperature recovery times after power outage and door openings.
- Better temperature uniformity

INCREASED SECURITY

Hold a set temperature inside the chamber much longer than air-jacketed units in the event of power failure.

COMPLETE CONTAMINATION CONTROL METHODS


- ULPA filter
- 0.2 micron inline filter
- ISOCIDE™ antimicrobial coating.
- 90°C Moist Heat Decon Cycle (water in the external chamber needs to be drained first)

CONVENIENCE

The unit is equipped with a fill-in port to place water and a drain valve to facilitate faster depletion of water when cleaning, before decon cycle or before transporting the equipment

EASY MONITORING

Water level can be check via the water level sensor.



OPTIONS AND ACCESSORIES

All options and accessories for standard CelCulture® CO₂ Incubators are also appropriate to use on Water-Jacketed CelCulture® CO₂ Incubators.

ORDERING INFORMATION

TC SENSOR MODEL WITH STAINLESS STEEL CHAMBER

Models		Description
230 V / 50-60 Hz	115 V / 50-60 Hz	
CCL-050A-8-WJ	CCL-050A-9-WJ	Celculture® Incubator, 50 L, TC Sensor, CO ₂ Control, Water-Jacketed
CCL-170A-8-WJ	CCL-170A-9-WJ	Celculture® Incubator, 170 L, TC Sensor, CO ₂ Control, ULPA, Water-Jacketed
CCL-240A-8-WJ	CCL-240A-9-WJ	Celculture® Incubator, 240 L, TC Sensor, CO ₂ Control, ULPA, Water-Jacketed

IR SENSOR MODEL WITH STAINLESS STEEL CHAMBER

Models		Description
230 V / 50-60Hz	115 V / 50-60 Hz	
CCL-050B-8-WJ	CCL-050B-9-WJ	Celculture® Incubator, 50 L, IR Sensor, CO ₂ Control, Water-Jacketed
CCL-170B-8-WJ	CCL-170B-9-WJ	Celculture® Incubator, 170 L, IR Sensor, CO ₂ Control, ULPA, Water-Jacketed
CCL-240B-8-WJ	CCL-240B-9-WJ	Celculture® Incubator, 240 L, IR Sensor, CO ₂ Control, ULPA, Water-Jacketed

5

SUPPRESSED O₂ MODEL WITH STAINLESS STEEL CHAMBER

Models		Description
230 V / 50-60 Hz	115 V / 50-60 Hz	
CCL-050T-8-WJ	CCL-050T-9-WJ	Celculture® Incubator, 50 L, IR Sensor, CO ₂ Control, O ₂ Control, Water-Jacketed
CCL-170T-8-WJ	CCL-170T-9-WJ	Celculture® Incubator, 170 L, IR Sensor, CO ₂ Control, O ₂ Control, ULPA, Water-Jacketed
CCL-240T-8-WJ	CCL-240T-9-WJ	Celculture® Incubator, 240 L, IR Sensor, CO ₂ Control, O ₂ Control, ULPA, Water-Jacketed

GENERAL SPECIFICATION

Models	CCL-50-_-WJ	CCL-170-_-WJ	CCL-240-_-WJ
Temperature			
Ambient Temp Range	18 to 34°C (64 to 93°F)		
Temperature Control Method	Direct Heat and Water Jacketed using PID microprocessor		
Temperature Range	Ambient +3°C to 60°C		
Temperature Uniformity	± 0.2°C*	± 0.2°C*	± 0.3°C*
Temperature Fluctuation	± 0.1°C	± 0.1°C	± 0.1°C
Temperature Recovery time** (after 1 min. door opening, 98% from initial value)	4 minutes	6 minutes	6 minutes
Start up time (25°C amb. to 37.0°C)	60 minutes	60 minutes	80 minutes
Power Off Temperature Drop Rate: 1 hour 10 hours	1.8°C 7.0°C	0.6°C 6.2°C	1.2°C 5.5°C
CO₂			
CO ₂ Control System	Microprocessor PID		
CO ₂ Range	0-20%		
CO ₂ Accuracy	± 0.1%		
CO ₂ Sensor	Infrared (IR) Sensor**/TC Sensor		
CO ₂ Recovery Time*** (after 1 min. door opening, 98% from initial value)	Standard Unit: 8 minutes Suppressed O ₂ model: 6 minutes	Standard Unit: 4 minutes Suppressed O ₂ model: 5 minutes	Standard Unit: 5 minutes Suppressed O ₂ model: 5 minutes
O₂ for Supressed O₂ Model			
O ₂ Control System	Microprocessor PID		
O ₂ Range	1% - 20.7%		
O ₂ Accuracy	± 0.1%		
O ₂ Sensor	Galvanic Cell Type		
O ₂ Recovery Time*** (after 1 min. door opening, 98% from initial value)	At 1.0% O ₂ volume: 10 minutes At 5.0% O ₂ volume: 6 minutes	At 1.0% O ₂ volume: 20 minutes At 5.0% O ₂ volume: 10 minutes	At 1.0% O ₂ volume: 24 minutes At 5.0% O ₂ volume: 12 minutes

6

Humidity			
Humidification Method	Humidity pan		
Humidity range*****	85% - 93%		
Humidity range (Suppressed O ₂ control)*****	85% - 91%		
Physical Parameters			
Interior Volume	50 L (1.8 cu. ft.)	170 L (6.0 cu. ft.)	240 L (8.5 cu. ft.)
Internal Dimensions (W x D x H)	345 mm x 375 mm x 390 mm (13.6" x 14.8" x 15.4")	505 mm x 530 mm x 635 mm (19.9" x 20.9" x 25.0")	595 mm x 620 mm x 635 mm (23.4" x 24.4" x 25.0")
External Dimension (W x D x H)	500 mm x 500 mm x 705 mm (19.7" x 19.7" x 27.8")	660 mm x 670 mm x 980 mm (26.0" x 26.4" x 38.6")	750 mm x 765 mm x 980 mm (29.5" x 30.1" x 38.6")
Water Jacket Volume	10 L	20 L	45 L
Net Weight	75 kg (165 lbs) (no water)	102 kg (225 lbs) (no water)	170 kg (374 lbs) (no water)
Shipping Weight	90 kg (198 lbs)	118 kg (260 lbs)	185 kg (407 lbs)
Shipping Dimension (W x D x H)	660 mm x 690 mm x 980 mm (26.0" x 27.2" x 38.6")	850 mm x 750 mm x 1240 mm (33.5" x 29.5" x 48.8")	860 mm x 860 mm x 1240 mm (33.8" x 33.8" x 48.8")
Number of Shelves	2	4	
Maximum No. of Shelves	4	7	
Shelves Area (W x D)	310 mm x 310 mm (12.2" x 12.2")	470 mm x 470 mm (18.5" x 18.5")	550 mm x 550 mm (21.7" x 21.7")
Max. Load per Shelf	4 kg/shelf (8.8 lbs/shelf)	11 kg/shelf (24.3 lbs/shelf)	15 kg/shelf (33.1 lbs/shelf)
Available Electrical Configuration	220 - 240 VAC, 50 / 60 Hz, 1 Φ , 3.4 A 110 - 130 VAC, 50 / 60 Hz, 1 Φ , 7.0 A		
Airflow	6-8 cfm		
Interior Material	Stainless steel, type 304		
Nominal Power 37°C	28 Watts	50 Watts	80 Watts
Power/FLA current	495 Watts, 2.8 A, 230 V 495 Watts, 5.6 A, 115 V	640 Watts, 2.8 A, 230 V 640 Watts, 5.6 A, 115 V	850 Watts, 3.7 A, 230 V 850 Watts, 7.4 A, 115 V
Contamination Control			
Contamination Control Methods	1) Main body is electrogalvanized steel with ISOCIDE™ antimicrobial coating; 2) 0.2 micron in-line filter for gas inputs; 3) ULPA filter***** 4) 90°C Moist Heat Decon Cycle (water in the external chamber needs to be drained first)		

* Data recorded under optimum factory setting conditions
 ** For temperature not exceeding 37°C
 *** For CO₂ not exceeding 5.2%. Recovery time for TC sensor is longer
 **** For O₂ not exceeding 5.2%
 ***** Esco does not guarantee condensation free chamber at higher humidity level.
 ***** Not available on 50 L size model

ESCO

WORLD CLASS. WORLDWIDE.



*CelCulture® CO₂ Incubator
Model CCL-170-_-*

CelCulture®

CO₂ Incubator with Integrated Cooling System

Solution for Highly Specialized Applications





CelCulture®

CO₂ Incubator with Integrated Cooling System

INTRODUCTION

Esco CelCulture CO₂ Incubator with Integrated Cooling System provides solution for highly specialized application.

The integrated cooling system allows studies of samples that requires temperature at/or below ambient temperature.

KEY FEATURES

WIDER TEMPERATURE RANGE

Temperature range of 12°C below ambient to 60°C above ambient means wider range of applications.

HIGHLY EFFICIENT, ENVIRONMENTALLY FRIENDLY PELTIER COOLING SYSTEM


This provides precise heating and cooling inside the chamber making sure that your samples are safe from temperature changes.

COMPLETE CONTAMINATION CONTROL METHODS

- 90 °C validated moist heat decontamination cycle
- ULPA filter
- ISOCIDE® anti-microbial coating
- 0.2 micron in-line filter

APPLICATIONS

- Stem Cell Research
- Mammalian Research
- Tissues Re-generation
- Food Research
- Regenerative Medicine
- Marine Biology
- Fish Cell Research
- Amphibian cell Research
- Insect Cell Research



CelCulture®

CO₂ Incubators with Integrated Cooling System

HEATING AND COOLING FUNCTION

Users can use it as a conventional CO₂ incubator using heating mode or switch to cooling mode.



```
SYSTEM SETPOINTS
TEMP=37.0°C  O2=5.0%
CO2= 5.0%  HEATING
SYSTEM OK
```



```
SYSTEM SETPOINTS
TEMP=12.0°C  O2=5.0%
CO2= 5.0%  COOLING
SYSTEM OK
```

OPTIONS AND ACCESSORIES

All options and accessories for standard CelCulture CO₂ incubators are also appropriate to use on CelCulture CO₂ Incubator with Integrated Cooling Option.

ORDERING INFORMATION

IR SENSOR MODEL WITH INTEGRATED COOLING SYSTEM

Models	Description
CCL-170-B-8-P	CelCulture Incubator, 170L, IR sensor, CO ₂ control, Moist Heat Decon, Peltier System, 230VAC, 50/60 Hz
CCL-170-B-9-P	CelCulture Incubator, 170L, IR sensor, CO ₂ control, Moist Heat Decon, Peltier System, 115VAC, 50/60 Hz

SUPPRESSED O₂ MODEL WITH INTEGRATED COOLING SYSTEM

Models	Description
CCL-170-T-8-P	CelCulture Incubator, 170L, IR sensor, CO ₂ control, O ₂ control, Moist Heat Decon, Peltier System, 230VAC, 50/60 Hz
CCL-170-T-9-P	CelCulture Incubator, 170L, IR sensor, CO ₂ control, O ₂ control, Moist Heat Decon, Peltier System, 115VAC, 50/60 Hz

GENERAL SPECIFICATION

Models	CCL-170B-8-P CCL-170B-9-P
TEMPERATURE (Heating Mode)	
Ambient Temp Range	18 to 34°C (64 to 93 °F)
Temperature Control Method	Direct Heat and Air Jacketed with Peltier Cooling System using PID microprocessor
Temperature Range, °C	12°C below ambient to 60°C above ambient
Temperature Uniformity, °C	± 0.2*
Temperature Fluctuation, °C	± 0.1
Temperature Recovery time** (after 1 min. door opening, 98% from initial value), minutes	6
TEMPERATURE (Cooling Mode)	
Pull Down Time to 8°C below Ambient (at 25°C Ambient Temperature), minutes	42
Pull Down Time to 12°C below Ambient (at 25°C Ambient Temperature), minutes	90
Temperature Uniformity at 8°C below Ambient, °C	± 0.4
Temperature Recovery Time*** (after 1 min. door opening, 98% from initial value), minutes	32
Power Off Temperature Increase Rate 1 hour, °C 10 hours, °C	1.9 10
CO ₂ (Heating Mode)	
CO ₂ Control System	Microprocessor PID
CO ₂ Range, % CO ₂	0-20
CO ₂ Accuracy, % CO ₂	± 0.1
CO ₂ Sensor	Infrared (IR) Sensor**
CO ₂ Recovery Time*** (after 1 min. door opening, 98% from initial value) (Heating Mode)	Standard Unit: 4 minutes Suppressed O ₂ model: 5 minutes
CO ₂ (Cooling Mode)	
Start-up Time at 5% CO ₂ , minutes	14
Recovery Time at 5% CO ₂ , minutes	9
CO ₂ Fluctuation, %CO ₂	± 0.3
FOR SUPPRESSED O ₂ MODEL (Heating Mode)	
O ₂ Control System	Microprocessor PID
O ₂ Range, % O ₂	1-20.7
O ₂ Accuracy, % O ₂	± 0.1
O ₂ Sensor	Galvanic Cell Type
O ₂ Recovery Time**** (after 1 min. door opening, 98% from initial value), minutes	At 1.0% O ₂ volume: 20 At 5.0% O ₂ volume: 10

6

FOR SUPPRESSED O ₂ MODEL (Cooling Mode)	
Start-up Time at 5% O ₂ , minutes	14
O ₂ Recovery Time at 5% O ₂ **** (after 1 min. door opening, 98% from initial value), minutes	12
O ₂ Fluctuation, %O ₂	± 0.2
Humidity (Heating Mode)	
Humidification Method	Humidity pan
Humidity range	85% - 93%
Humidity range (Suppressed O ₂ control)	85% - 91%
Physical Parameters	
Interior Volume	170l (5.7 cu. Ft.)
Internal Dimensions (W x D x H)	505 mm x 530 mm x 635 mm (19.9" x 20.9" x 25.0")
External Dimension (W x D x H)	660 mm x 672 mm x 900 mm (26.0" x 26.46" x 35.4")
Net Weight	105 kg (231 lbs)
Shipping Weight	131.5 kg (289 lbs)
Shipping Dimension (W x D x H)	820 mm x 850 mm x 930 mm 32.28" x 33.47" x 36.62"
Number of Shelves	4
Maximum No. of Shelves	7
Shelves Area (W x D)	470 mm x 470 mm (18.5" x 18.5")
Max. Load per Shelf	11 kg/shelf (24.3 lbs/shelf)
Available Electrical Configuration	220 - 240 VAC, 50 / 60 Hz, 1Φ, 3.4 A 110 - 130 VAC, 50 / 60 Hz, 1Φ, 7.0 A
Airflow	6-8 cfm
Interior Material	Stainless steel, type 304
Maximum Power Consumption	800 Watts
Power Consumption 37°C	80 W
CONTAMINATION CONTROL	
Contamination Control Methods	1) Main body is electrogalvanized steel with ISOCIDE® antimicrobial coating; 2) Moist 90°C OVERNIGHT decon. cycle (HPA validated); 3) 0.2 micron in-line filter for gas inputs; 4) ULPA filter

* Data recorded under optimum factory setting conditions
 ** For temperature not exceeding 37°C
 *** For CO₂ not exceeding 5.2%.
 **** For O₂ not exceeding 5.2%



WORLD CLASS. WORLDWIDE.

CelCulture®

CO₂ Incubator with Stainless Steel Exterior Cabinet

- Corrosion Resistant Surface
- Meets Pharmaceutical and Clinical Laboratory Requirements



ORDERING INFORMATION

TC SENSOR MODEL WITH STAINLESS STEEL EXTERIOR CABINET

Models	Description
CCL-050A-8-SS	CelCulture® Incubator, 50 L, TC sensor, CO ₂ Control, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz
CCL-170A-8-SS	CelCulture® Incubator, 170 L, TC Sensor, CO ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz
CCL-240A-8-SS	CelCulture® Incubator, 240 L, TC Sensor CO ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz

IR SENSOR MODEL WITH STAINLESS STEEL EXTERIOR CABINET

Models	Description
CCL-050B-8-SS	CelCulture® Incubator, 50 L, IR sensor, CO ₂ Control, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz
CCL-170B-8-SS	CelCulture® Incubator, 170 L, IR Sensor, CO ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz
CCL-240B-8-SS	CelCulture® Incubator, 240 L, IR Sensor CO ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz

SUPPRESSED O₂ MODEL WITH STAINLESS STEEL EXTERIOR CABINET

Models	Description
CCL-050T-8-SS	CelCulture® Incubator, 50 L, IR sensor, CO ₂ & O ₂ Control, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz
CCL-170B-8-SS	CelCulture® Incubator, 170 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz
CCL-240B-8-SS	CelCulture® Incubator, 240 L, IR Sensor, CO ₂ & O ₂ Control, ULPA, Moist Heat Decon, SS Cabinet, 230 VAC, 50/60 Hz

GENERAL SPECIFICATIONS

CELCULTURE® CO₂ INCUBATORS

	CCL-050_ _ -SS	CCL-170_ _ -SS	CCL-240_ _ -SS
Temperature			
Temperature Control Method	Direct heat & air jacket using Microprocessor PID		
Temperature Range	Amb. +3°C to 60°C		
Temperature Uniformity	<± 0.2°C*	<± 0.2°C*	<± 0.3°C*
Temperature Accuracy	<± 0.1°C		
Recovery Time** (after 1 min. door opening, 98% from initial value)	4 minutes	6 minutes	6 minutes
Ambient Temperature Range	18 to 34°C (64 to 93°F)		
CO₂			
CO ₂ Control System	Microprocessor PID		
CO ₂ Range	0-20%		
CO ₂ Accuracy	± 0.1%		
CO ₂ Sensor	Infrared (IR) Sensor*** / TC Sensor		
CO ₂ Recovery Time**** (after 1 min. door opening, 98% from initial value)	Standard Unit: 8 minutes Suppressed O ₂ model: 8 minutes	Standard Unit: 4 minutes Suppressed O ₂ model: 5 minutes	Standard Unit: 5 minutes Suppressed O ₂ model: 5 minutes
O₂ Specs (For Suppressed O₂ Model)			
O ₂ Control System	Microprocessor PID		
O ₂ Range	1-20.7%		
O ₂ Accuracy	± 0.1%		
O ₂ Sensor	Galvanic Cell Type		
O ₂ Recovery Time (after 1 minute door opening)	at 1.0% O ₂ by volume: 10 minutes	at 1.0% O ₂ by volume: 20 minutes	at 1.0% O ₂ by volume: 24 minutes
	at 5.0% O ₂ by volume: 6 minutes	at 5.0% O ₂ by volume: 10 minutes	at 5.0% O ₂ by volume: 12 minutes
Humidity			
Humidification Method	Humidity pan		
Humidity Range*****	85% - 93% RH		
Humidity Range (Suppressed O ₂ Control)*****	85% - 91%		
Physical Construction			
Interior Volume	50 L (1.8 cu. ft.)	170 L (6.0 cu. ft.)	240 L (8.5 cu. ft.)
External Dimensions (W x D x H)	500 mm x 500 mm x 655 mm (19.7" x 19.7" x 25.8")	660 mm x 660 mm x 900 mm (26.0" x 26.0" x 35.4")	750 mm x 665 mm x 900 mm (29.5" x 26.2" x 35.4")
Internal Dimensions (W x D x H)	345 mm x 375 mm x 390 mm (13.6" x 14.8" x 15.4")	505 mm x 530 mm x 635 mm (19.9" x 20.9" x 25.0")	595 mm x 620 mm x 635 mm (23.4" x 24.4" x 25.0")
Shipping Weight	70 kg (154.3 lbs)	120 kg (264.6 lbs)	155 kg (341.7 lbs)
Shipping Dimensions (W x D x H)	660 mm x 660 mm x 890 mm (26.0" x 26.0" x 35.0")	850 mm x 720 mm x 1150 mm (33.5" x 28.3" x 45.3")	860 mm x 830 mm x 1110 mm (33.9" x 32.7" x 43.7")
Number of Shelves	3	4	4
Maximum No. of Shelves	4	7	7
Shelves Area (W x D)	310 mm x 310 mm (12.2" x 12.2")	470 mm x 470 mm (18.5" x 18.5")	550 mm x 550 mm (21.7" x 21.7")
Max. Load per Shelf	4 kg/shelf (8.8 lbs/shelf)	11 kg/shelf (24.3 lbs/shelf)	15 kg/shelf (33.1 lbs/shelf)
Available Electrical Configuration	220 - 240 VAC, 50/60 Hz, 1ϕ, 3.4 A		
	110 - 130 VAC, 50/60 Hz, 1ϕ, 7.0 A		
Maximum Power Consumption	372 Watts	800 Watts	1100 Watts
Power Consumption 37°C	37 Watts	80 Watts	110 Watts
Interior Material	Stainless steel, type 304		
Contamination Control			
Contamination Control Methods	1) Moist 90°C OVERNIGHT decon cycle (HPA validated) 2) 0.2 micron in-line filter for gas inputs 3) ULPA filter*****		

* Data recorded under optimum factory setting conditions

** For temperature not exceeding 37°C

*** For CO₂ not exceeding 5.2%. Recovery time for TC sensor is longer

**** For O₂ not exceeding 5.2%

***** Esco does not guarantee condensation free chamber at higher humidity level.

***** Not available on 50 L size model



WORLD CLASS. WORLDWIDE.

Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777

Tel +65 6542 0833 • Fax +65 6542 6920 • mail@escoglobal.com

www.escoglobal.com

