



## For determination of Phosphorus

*Sample preparation mineralizer D 8P for determination of Phosphorus!*

**Mineralizer BEGER D 8P is designed to accelerate the determination of Phosphorus in plant materials, feed and fodder raw materials.**



**Mineralizer BEGER D 8P is designed to accelerate the determination of Phosphorus in plant materials, feed and fodder raw materials.**

**Field of usage – industrial laboratories of the factories that store and processing grain and cereal products, quality control laboratories of food products and research institutions.**

- ✓ **The newest device has no analogues that provide the mineralization process**
- ✓ **Analysis for up to 8 samples of 50 ml simultaneously**
- ✓ **Saving of electricity, working time and reagents**
- ✓ **Temperature control**

## **Specifications:**

<b>Number of samples</b>	<b>8</b>
<b>Dimensions (W x D x H)</b>	<b>390x400x280 mm</b>
<b>Control</b>	<b>Microprocessor PID controller</b>
<b>Power</b>	<b>800 W</b>
<b>Heating range</b>	<b>25 to 400 °C</b>
<b>Flasks volume</b>	<b>50 ml</b>
<b>Weight</b>	<b>7 kg</b>
<b>Cat. N°</b>	<b>514 000 00 01</b>



**Stegne 27, Industrijska cona Stegne,  
1000 Ljubljana, Slovenia**

**[www.beger.si](http://www.beger.si)**

**e-mail: [info@beger.si](mailto:info@beger.si)**

**skype: [beger\\_d.o.o](https://www.skype.com/people/beger_d.o.o)**

**Find us on:**

